

Assembly Instructions

AML25059 - Raised 2IN Toilet Seat

AML25061 - Raised 4IN Toilet Seat

Check you have all the parts:

| | |
|----|----------------|
| 1x | one-piece seat |
|----|----------------|

Tools required:

| |
|-----|
| N/A |
|-----|

Instructions

1. To fasten to toilet seat, loosen the butterfly screws on the 3 mounts
2. Place seat on toilet and tighten to suit the toilet seat size

Guidance for Use

Maximum Weight Capacity: 190kg