

Q6 EDGE 3 *Stretto*™



Q6 EDGE 3 STRETTO™
shown with iLevel® power seat elevation,
Q-Logic 3 environmental controls, and
Stealth® TRU-Comfort 2 cushion.



The Q6 Edge 3 Stretto™ features a narrow width for outstanding manoeuvrability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 12.5" (31.8cm) or 14" (35.6cm) drive wheels and optional 12" (30.5cm) of iLevel® power adjustable seat height.



Narrowest
complex
rehab base
on the market


QUANTUM®
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	14" (35.6cm) or 12.5" (31.8cm) pneumatic or flat free tyres
CASTOR WHEELS	Front: 5" (12.7cm) flat free Rear: 6" (15.2cm) flat free
MAXIMUM SPEED ¹	Up to 10 km/h Up to 5.6 km/h w/iLevel activated
GROUND CLEARANCE	2.9" (7.4cm)
TURNING RADIUS ²	21" (53.3cm)
OVERALL LENGTH ²	32.5" (82.5cm) without front riggings
BASE WIDTH	14" Drive wheels: 21.75" (55.2cm) 12.5" Drive wheels: 20.47" (52cm)
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 12-20" (30.5cm-51cm) Depth: 12-20" (30.5cm-51cm)
TRU-Balance® 3 Power Tilt	Width: 12-20" (30.5cm-51cm) Depth: 12-20" (30.5cm-51cm)
TRU-Balance® 3 Tilt and iLevel®	Width: 12-20" (30.5cm-51cm) Depth: 12-20" (30.5cm-51cm)
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	16-18" (40.6cm - 45.7cm)
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5" (44.5cm - 49.5cm)
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5-19.5" (44.5cm- 49.5cm)
BATTERY SIZE ³	55 AH (2)
BATTERY WEIGHT ⁶	55 AH: 18.1 kg each
AVAILABLE ELECTRONICS ⁷	120A Q-Logic 3 EX
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6®
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT WITHOUT SEATING SYSTEM	57.3 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 0.3 cm with 45.7 cm x 45.7 cm seat, fixed towers, 136 kg weight capacity and solid tyres (add 0.95 cm with pneumatic tyres). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) Battery weight may vary +/- 1.4 kg based on manufacturer.
 7) All Quantum Rehab® units can be ordered with specialty controls.

Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 30.5 cm of seat elevation at walking speed (up to 5.6 km/h)
- Fender lights and USB charger included
- Base made in America of imported and domestic parts
- Assembled to specification in New Zealand

Options

- iLevel® option allows 30.5 cm of seat elevation while at walking speed (5.6 km/h)
- TRU-Balance® 3 Seating
- 30.5 cm x 30.5 cm seat sizes with iLevel
- Back in Black shrouds (Green, blue and raspberry shrouds available for indent orders)

